

SCHWEIZER Magnifier MiniMag

Operating Instructions Pocket Magnifier with LED light

CE

1 Intended use

The product is a magnifier. Use the product to enlarge objects.

The product is intended for indoor use only. Do not use it outdoors. Contact with moisture must be avoided under all circumstances.

If you use the product for purposes other than those described, the product may be damaged.

Improper use can result in short circuits, fires, or other hazards.

The product complies with the statutory national and European requirements.

For safety and approval purposes, you must not rebuild and/or modify the product. Read the operating instructions carefully and store them in a safe place. Make this product available to third parties only together with the operating instructions.

2 Delivery contents

- Product
- Operating instructions

3 Description of symbols

The following symbols are on the product/appliance or are used in the text:



The symbol warns of hazards that can lead to personal injury.



This product is constructed according to protection class III.



This product must only be used in dry, enclosed indoor areas. It must not become damp or wet.

4 Safety instructions



Read the operating instructions carefully and especially observe the safety information. If you do not follow the safety instructions and information on proper handling, we assume no liability for any resulting personal injury or damage to property. Such cases will invalidate the warranty/guarantee.

4.1 General

The product is not a toy. Keep it out of the reach of children and pets.

- Do not leave packaging material lying around carelessly. This may become dangerous playing material for children.
- Handle the product carefully. Jolts, impacts or a fall even from a low height can damage the product.

4.2 Risk of fire

- If the lens is exposed to sunlight, the magnifier lens acts as a burning lens. The heat can ignite combustible materials and result in fire. **Always cover the lens if not in use.**

4.3 Eye injury

- Never use the magnifier to stare directly into sunlight or other strong light sources. Do not use the magnifier near them. This poses a risk of loss of eyesight!

4.4 LED light

- Do not look directly into the LED light!
- Do not look into the beam directly or with optical instruments!

4.5 Operating environment

- Do not place the product under any mechanical stress.
- Protect the appliance from extreme temperatures, strong jolts, flammable gases, steam and solvents.
- Protect the product from high humidity and moisture.
- Protect the product from direct sunlight.

4.6 Operation

- Consult an expert when in doubt about the operation, safety or connection of the product. If it is no longer possible to operate the product safely, take it out of operation and protect it from any accidental use. **DO NOT** attempt to repair the product yourself. Safe operation can no longer be guaranteed if the product:
 - is visibly damaged,
 - is no longer working properly,
 - has been stored for extended periods in poor ambient conditions or
 - has been subjected to any serious transport-related stresses.

4.7 (Rechargeable) Batteries

- Correct polarity must be observed while inserting (rechargeable) batteries.
- The (rechargeable) batteries should be removed from the device if it is not used for a long period of time to avoid damage through leaking. Leaking or damaged (re-chargeable) batteries might cause acid burns when in contact with skin, therefore use suitable protective gloves to handle corrupted (rechargeable) batteries.
- (Rechargeable) batteries must be kept out of reach of children. Do not leave (re-chargeable) batteries lying around, as there is risk, that children or pets swallow them.
- All (rechargeable) batteries should be replaced at the same time. Mixing old and new (re-chargeable) batteries in the device can lead to (rechargeable) battery leakage and device damage.
- (Rechargeable) batteries must not be dismantled, short-circuited or thrown into fire. Never recharge non-rechargeable batteries. There is a risk of explosion!