

eMag 43 HD

Optimized image, operation and design



Technical data

eMag 43 HD

The portable video magnifier in optimized HD quality with folding handle can be used like a traditional magnifying glass or can also be positioned directly on the object.

- Item No. 544002
- Magnification 1.6x to 14x continuous zoom
- Viewing options Full colour (with quick change feature), black on white, white on black, yellow on blue, yellow on black + 15 additional contrast modes
- Weight 228 g
- Size Size (WxDxH) approx. 126x75x32 mm (closed)
- Operating time Approx. 2 hours
- Battery 4 rechargeable NiMH batteries AAA (approx. 2 hours charge time)
- Display 10,9 cm (4.3") LCD display 16:9
- Features
 - Snapshot function, 80 images can be saved and exported via USB
 - Reading line and text masks
 - LED battery indicator
 - LED illumination can be deactivated
 - Energy saving mode after 30 minutes
- Accessories
 - Case with belt loop and carrying strap
 - Wrist strap
 - AC power adapter
 - USB cable
 - Manual
- Extras Item 544042: Writing and reading stand



A. SCHWEIZER GmbH
 Hans-Böckler-Str. 7, 91301 Forchheim, Germany
 Phone: +49-9191-721056, Fax: +49-9191-721073
 info@schweizer-optik.de, www.schweizer-optik.com

23-041/313118/02