

eMag 43

The classic portable video magnifier



Technical data

eMag 43

The portable video magnifier with folding handle has a centrally positioned camera and can be used like a traditional magnifying glass or can also be positioned directly on the object.

- Item No. 543002
- Magnification 2x to 14x in 3 steps
- Viewing options Full colour, black on white, white on black, yellow on blue, yellow on black
- Weight 218 g
- Size Size (WxDxH) approx. 126x75x32 mm (closed)
- Operating time Approx. 2 hours
- Battery 4 rechargeable NiMH batteries AAA (approx. 2 hours charge time)
- Display 10.9 cm (4.3") LCD display 16:9
- Features
 - Snapshot function, a maximum of 15 snapshots can be saved
 - LED battery indicator
 - LED illumination can be deactivated
 - Energy saving mode after 30 minutes
- Accessories
 - Case with belt loop and carrying strap
 - Wrist strap
 - AC power adapter
 - Manual
- Extras Item 535042: Writing and reading stand



A. SCHWEIZER GmbH
Hans-Böckler-Str. 7, 91301 Forchheim, Germany
Phone: +49-9191-721056, Fax: +49-9191-721073
info@schweizer-optik.de, www.schweizer-optik.com

28-041/313288/00